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**Effectiveness of peripheral defocus spectacle lenses**  
(Ma et al., 2022)<sup>5</sup>

**D.I.M.S. Technology: 2-year RCT**  
(Lam et al., 2020)<sup>6</sup>

**Impact of D.I.M.S. Technology on visual function**  
(Lam et al., 2020)<sup>7</sup>

**Impact of D.I.M.S. Technology on relative peripheral refraction** (Zhang et al., 2020)<sup>8</sup>

**D.I.M.S. Technology: 3-year follow-up clinical study**  
(Lam et al., 2022)<sup>9</sup>

**D.I.M.S. Technology: 6-year follow-up clinical study**  
(Lam et al., 2023)<sup>10</sup>

**Performance of MiYOSMART spectacle lenses in European children** (McCullough et al., 2023)<sup>11</sup>

**Adaptation and acceptance study of D.I.M.S. Technology** (Lu et al., 2020)<sup>12</sup>

**Effectiveness of D.I.M.S. Technology in combination with atropine** (Nucci et al., 2023)<sup>13</sup>

**Effectiveness of D.I.M.S. Technology in diverse clinical circumstances**  
(Liu et al., 2023)<sup>14</sup>

**Effects of myopic defocus in chicks** (Tse et al., 2011)<sup>15</sup>

**Effects of myopic defocus in infant monkeys**  
(Arumugam et al., 2014)<sup>16</sup>

**Simultaneous myopic and hyperopic defocus in chicks**  
(Tse et al., 2007)<sup>17</sup>

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**Effectiveness of peripheral defocus spectacle lenses**  
(Ma et al., 2022)<sup>5</sup>

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**D.I.M.S. Technology:**  
**2-year RCT**  
(Lam et al., 2020)<sup>6</sup>

**Impact of D.I.M.S. Technology on visual function**  
(Lam et al., 2020)<sup>7</sup>

**Impact of D.I.M.S. Technology on relative peripheral refraction**  
(Zhang et al., 2020)<sup>8</sup>

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**D.I.M.S. Technology:**  
**3-year follow-up clinical study**  
(Lam et al., 2022)<sup>9</sup>

**D.I.M.S. Technology:**  
**6-year follow-up clinical study**  
(Lam et al., 2023)<sup>10</sup>

**Performance of MiYOSMART spectacle lenses in European children: 12-month outcomes** (McCullough et al., 2023)<sup>11</sup>

**Adaptation and acceptance study of D.I.M.S. Technology**  
(Lu et al., 2020)<sup>12</sup>

**Effectiveness of D.I.M.S. Technology in combination with atropine**  
(Nucci et al., 2023)<sup>13</sup>

**Effectiveness of D.I.M.S. Technology in diverse clinical circumstances**  
(Liu et al., 2023)<sup>14</sup>

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**Effects of myopic defocus in chicks**  
(Tse et al., 2011)<sup>15</sup>

**Effects of myopic defocus in infant monkeys**  
(Arumugam et al., 2014)<sup>16</sup>

**Simultaneous myopic and hyperopic defocus in chicks**  
(Tse et al., 2007)<sup>17</sup>

### Referenzen:

1. Brennan NA, et al. Efficacy in myopia control. Prog Retin Eye Res. 2021;83:100923.
2. Tricard D, et al. Progression of myopia in children and teenagers: a nationwide longitudinal study. Br J Ophthalmol. 2022;106(8):1104-9.
3. Morais F, et al. Levels of Evidence: What Should Ophthalmologists Know? Rev Bras Oftalmol. 2019;78(6):413-7.
4. Lawrenson JG, et al. Interventions for myopia control in children: a living systematic review and network meta-analysis. Cochrane Database Syst Rev. 2023;2(2):CD014758.
5. Ma JX, et al. Effectiveness of peripheral defocus spectacle lenses in myopia control: a Meta-analysis and systematic review. Int J Ophthalmol. 2022;15(10):1699-706.
6. Lam CSY, et al. Defocus Incorporated Multiple Segments (DIMS) spectacle lenses slow myopia progression: a 2-year randomised clinical trial. Br J Ophthalmol. 2020;104(3):363-8.
7. Lam CSY, et al. Effect of Defocus Incorporated Multiple Segments Spectacle Lens Wear on Visual Function in Myopic Chinese Children. Transl Vis Sci Technol. 2020;9(9):11.
8. Zhang HY, et al. Defocus Incorporated Multiple Segments Spectacle Lenses Changed the Relative Peripheral Refraction: A 2-Year Randomized Clinical Trial. Invest Ophthalmol Vis Sci. 2020;61(5):53.
9. Lam CS, et al. Myopia control effect of defocus incorporated multiple segments (DIMS) spectacle lens in Chinese children: results of a 3-year follow-up study. Br J Ophthalmol. 2022;106(8):1110-4.
10. Lam CSY, et al. Long-term myopia control effect and safety in children wearing DIMS spectacle lenses for 6 years. Sci Rep. 2023;13(1):5475.
11. McCullough S, et al. Presented at: The Association for Research in Vision and Ophthalmology, 23-27 April 2023, New Orleans, LA, USA (Abstract).
12. Lu Y, et al. The Adaptation and Acceptance of Defocus Incorporated Multiple Segment Lens for Chinese Children. Am J Ophthalmol. 2020;211:207-16.
13. Nucci P, et al. A comparison of myopia control in European children and adolescents with defocus incorporated multiple segments (DIMS) spectacles, atropine, and combined DIMS/atropine. PLoS One. 2023;18(2):e0281816.
14. Liu J, et al. The Efficacy of Defocus Incorporated Multiple Segments Lenses in Slowing Myopia Progression: Results from Diverse Clinical Circumstances. Ophthalmology. 2023;130(5):542-50.
15. Tse DY, To CH. Graded competing regional myopic and hyperopic defocus produce summated emmetropization set points in chick. Invest Ophthalmol Vis Sci. 2011;52(11):8056-62.
16. Arumugam B, et al. The effects of simultaneous dual focus lenses on refractive development in infant monkeys. Invest Ophthalmol Vis Sci. 2014;55(11):7423-32.
17. Tse DY, et al. Simultaneous defocus integration during refractive development. Invest Ophthalmol Vis Sci. 2007;48(12):5352-9.

**Mehr Informationen  
zu MiYOSMART**

**HOYA**  
FOR THE VISIONARIES